



Item Number M1148

Item Methylcellulose, 400 cP, FCC

CAS Number 9004-67-5

Molecular Formula

Molecular Weight

Assay [Methoxyl Groups (-OCH₃); Dried Basis]..... 27.5-31.5%

Viscosity of a 2 g in 100 ml Solution..... 300-560 cP

MAXIMUM LIMITS

Lead..... 3 mg/kg

Loss on Drying..... 5.0%

Residue on Ignition..... 1.5%

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: ***East Coast Facility:***

14422 S. San Pedro St. 755 Jersey Ave.

Gardena, CA 90248 New Brunswick, NJ 08901

(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843