



Item Number M1198

Item Mineral Oil, White, Heavy, FCC

CAS Number 8012-95-1

Molecular Formula

Molecular Weight

Carbon number at 5% distillation point..... Min. 28

Specific Gravity.....0.845-0.905

Ultraviolet Absorbance (Polynuclear Hydrocarbons)..... To pass test

Viscosity @ 100°C..... Min. 11 cSt

MAXIMUM LIMITS

Lead..... 1 mg/kg

Readily Carbonizable Substances..... To pass test

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: East Coast Facility:

14422 S. San Pedro St. 755 Jersey Ave.

Gardena, CA 90248 New Brunswick, NJ 08901

(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843