



Item Number M1243
Item Methyl Benzoate
CAS Number 93-58-3
Molecular Formula C₈H₈O₂
Molecular Weight 136.14

Assay (C₈H₈O₂)..... Min. 98%
Specific Gravity @ 25°C..... 1.082-1.088
Refractive Index @ 20°C..... 1.514-1.518
Solubility in Alcohol..... To pass test
MAXIMUM LIMITS
Acid Value..... 1.0
Chlorinated Compounds..... To pass test

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: East Coast Facility:

14422 S. San Pedro St. 755 Jersey Ave.

Gardena, CA 90248 New Brunswick, NJ 08901

(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843