



Item Number M1277
Item Methylparaben Sodium, NF
CAS Number 5026-62-0
Molecular Formula $C_8H_7NaO_3$
Molecular Weight 174.14

Assay (Anhydrous Basis).....98.5-101.5%
Completeness of Solution..... To pass test
pH of a 1 in 1000 Solution..... 9.5-10.5

MAXIMUM LIMITS

Water..... 5.0%
Chloride.....0.035%
Sulfate.....0.12%
Organic Volatile Impurities..... To pass test

Spectrum Chemicals & Laboratory Products

<i>Corporate Headquarters:</i>	<i>East Coast Facility:</i>
14422 S. San Pedro St.	755 Jersey Ave.
Gardena, CA 90248	New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843	