



Item Number M3239
Item Mitotane, USP
CAS Number 53-19-0
Molecular Formula $C_{14}H_{10}Cl_4$
Molecular Weight 320.04
Assay (Dried Basis)..... 97.0-103.0%
Melting Range..... 75°-81°C
MAXIMUM LIMITS
Loss on Drying..... 0.5%
Residue on Ignition..... 0.5%
Residual Solvents..... To pass test

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: East Coast Facility:
14422 S. San Pedro St. 755 Jersey Ave.
Gardena, CA 90248 New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843