



Item Number ME170
Item Methyl Salicylate, NF
CAS Number 119-36-8
Molecular Formula $C_8H_8O_3$
Molecular Weight 152.15

Assay.....98.0-100.5%
Solubility in 70% Alcohol..... To pass test
Specific Gravity @ 25/25°C..... 1.180-1.185
Angular Rotation..... Optically inactive
Refractive Index @ 20°C..... 1.535-1.538

MAXIMUM LIMITS

Heavy Metals..... 20 µg/g
Residual solvents..... To pass test

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: East Coast Facility:

14422 S. San Pedro St. 755 Jersey Ave.

Gardena, CA 90248 New Brunswick, NJ 08901

(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843