



Item Number ME175
Item Methyl Violet 2B
CAS Number 8004-87-3
Molecular Formula C₂₄H₂₈N₃Cl
Molecular Weight 393.96
Absorption Max. (in 50% Ethanol)..... 584 nm
Visual Transition Interval.....pH 0.1 (yellow) to pH 2.7 (violet)

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: 14422 S. San Pedro St.
Gardena, CA 90248
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843

East Coast Facility: 755 Jersey Ave.
New Brunswick, NJ 08901