



Item Number P1013
Item Pentane, Reagent
CAS Number 109-66-0
Molecular Formula C₅H₁₂
Molecular Weight 72.15

Assay (by GC)..... Min. 98%

MAXIMUM LIMITS

Color (APHA).....10

Residue after Evaporation..... 0.002%

Water..... 0.02%

Sulfur Compounds (as S)..... 0.005%

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: ***East Coast Facility:***

14422 S. San Pedro St. 755 Jersey Ave.

Gardena, CA 90248 New Brunswick, NJ 08901

(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843