



Item Number P1315
Item Potassium Hydroxide, Pellets, Reagent, ACS
CAS Number 1310-58-3
Molecular Formula KOH
Molecular Weight 56.11

Assay (KOH)..... Min. 85%

MAXIMUM LIMITS

Potassium Carbonate (K_2CO_3).....2.0%
Chloride (Cl)..... 0.01%
Nitrogen Compounds (as N)..... 0.001%
Phosphate (PO_4).....5 ppm
Sulfate (SO_4).....0.003%
Heavy Metals (as Ag).....0.001%
Iron (Fe).....0.001%
Nickel (Ni).....0.001%
Calcium (Ca)..... 0.005%
Magnesium (Mg).....0.002%
Sodium (Na)..... 0.05%

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: East Coast Facility:

14422 S. San Pedro St. 755 Jersey Ave.

Gardena, CA 90248 New Brunswick, NJ 08901

(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843