



Item Number P1346
Item Potassium Metabisulfite, Crystal, NF
CAS Number 16731-55-8
Molecular Formula $K_2S_2O_5$
Molecular Weight 222.33
Assay (as SO_2)..... 51.8-57.6%

MAXIMUM LIMITS

Iron..... 0.001%
Heavy Metals.....0.001%
Residual solvents..... To pass test

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: ***East Coast Facility:***
14422 S. San Pedro St. 755 Jersey Ave.
Gardena, CA 90248 New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843