



Item Number P1395

Item Potassium Hydrogen Sulfate, Fused, Powder,

CAS Number 7790-62-7

Molecular Formula

Molecular Weight

Acidity (as H_2SO_4)..... 37.5-38.6%

MAXIMUM LIMITS

Water (H_2O)..... 2.5%

Insoluble Matter..... 0.01%

Chloride (Cl)..... 0.002%

Phosphate (PO_4)..... 0.001%

Heavy Metals (as Pb)..... 0.001%

Iron (Fe).....0.002%

Calcium..... 0.002%

Magnesium..... 0.001%

Sodium (Na)..... 0.01%

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: East Coast Facility:

14422 S. San Pedro St. 755 Jersey Ave.

Gardena, CA 90248 New Brunswick, NJ 08901

(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843