



Item Number	P1402
Item	Sodium Potassium Tartrate, Tetrahydrate,
CAS Number	6381-59-5
Molecular Formula	$C_4H_4KNaO_6 \cdot 4H_2O$
Molecular Weight	282.22
	Assay ($KNaC_4H_4O_6 \cdot 4H_2O$)..... 99.0-102.0%
	Alkalinity..... To pass test
	Water..... 21.0-26.0%
	MAXIMUM LIMITS
	Lead..... 2 mg/kg

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: ***East Coast Facility:***

14422 S. San Pedro St. 755 Jersey Ave.

Gardena, CA 90248 New Brunswick, NJ 08901

(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843