



Item Number PA115

Item Pancreatin, Powder, USP

CAS Number 8049-47-6

Molecular Formula

Molecular Weight

Assay:

Amylase Activity..... Min. 25 U.S.P. units/mg

Lipase Activity..... Min. 2.0 U.S.P. units/mg

Protease Activity..... Min. 25 U.S.P. units/mg

MAXIMUM LIMITS

Loss on Drying..... 5.0%

Microbial Limits:

Salmonella..... To pass test

Escherichia coli..... To pass test

Fat..... 3.0%

Residual solvents..... To pass test

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: East Coast Facility:

14422 S. San Pedro St. 755 Jersey Ave.

Gardena, CA 90248 New Brunswick, NJ 08901

(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843