



Item Number PH125
Item Liquefied Phenol, USP
CAS Number 108-95-2
Molecular Formula C_6H_6O
Molecular Weight 94.11

Assay (C_6H_6O ; by Wt.)..... Min. 89.0%
Clarity of Solution and Reaction..... To pass test

MAXIMUM LIMITS

Distilling Range..... 182.5°C
Nonvolatile Residue..... 0.05%
Residual Solvents..... To pass test

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: ***East Coast Facility:***

14422 S. San Pedro St. 755 Jersey Ave.

Gardena, CA 90248 New Brunswick, NJ 08901

(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843