



Item Number PO170
Item Potassium Chloride, Granular, USP
CAS Number 7447-40-7
Molecular Formula KCl
Molecular Weight 74.55
Assay.....99.0-100.5%

MAXIMUM LIMITS

Acidity or Alkalinity..... To pass test
Loss on Drying..... 1.0%
Iodide or Bromide..... To pass test
Calcium and Magnesium..... To pass test
Sodium..... To pass test
Heavy Metals.....0.001%
Residual solvents..... To pass test

Spectrum Chemicals & Laboratory Products

<i>Corporate Headquarters:</i>	<i>East Coast Facility:</i>
14422 S. San Pedro St.	755 Jersey Ave.
Gardena, CA 90248	New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843	