



Item Number S1055
Item Silver Chloride, Granular, Reagent
CAS Number 7783-90-6
Molecular Formula AgCl
Molecular Weight 143.32
Assay (AgCl)..... Min. 99.0%

MAXIMUM LIMITS

Copper (Cu)..... 0.001%
Lead (Pb)..... 0.005%
Iron (Fe)..... 0.001%
Nitrate (NO₃)..... 0.001%

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: ***East Coast Facility:***
14422 S. San Pedro St. 755 Jersey Ave.
Gardena, CA 90248 New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843