



Item Number S1370
Item Sodium Oxalate, Powder, Reagent, ACS
CAS Number 62-76-0
Molecular Formula $C_2Na_2O_4$
Molecular Weight 134.00

Assay [$C_2Na_2O_4$]..... Min. 99.5%
Neutrality..... To pass test

MAXIMUM LIMITS

Insoluble Matter..... 0.005%
Loss on Drying..... 0.01%
Chloride (Cl)..... 0.002%
Ammonium (NH_4)..... 0.002%
Sulfate (SO_4)..... 0.002%
Heavy Metals (as Pb)..... 0.002%
Iron (Fe)..... 0.001%
Potassium (K)..... 0.005%
Substances Darkened by Hot Sulfuric Acid..... To pass test

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: East Coast Facility:

14422 S. San Pedro St. 755 Jersey Ave.

Gardena, CA 90248 New Brunswick, NJ 08901

(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843