



Item Number S1419
Item Sodium Pyrophosphate, Anhydrous, Powder,
CAS Number 7722-88-5
Molecular Formula $\text{Na}_4\text{P}_2\text{O}_7$
Molecular Weight 265.90
Assay ($\text{Na}_4\text{P}_2\text{O}_7$; Ignited Basis).....95.0-100.5%
MAXIMUM LIMITS
Arsenic (As)..... 3 mg/kg
Fluoride (F)..... 0.005%
Insoluble Matter..... 0.2%
Lead..... 4 mg/kg
Loss on Ignition.....0.5%

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: East Coast Facility:
14422 S. San Pedro St. 755 Jersey Ave.
Gardena, CA 90248 New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843