

Date Prepared: 03/01/2011

Reviewed On: 02/11/2011

1 Identification of the substance/mixture and of the company/undertaking

- **Product Identifier**
- **Product Name:** Dermatophyte Test Medium Base
- **Catalog Number:** 212330
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
BD Diagnostic Systems
7 Loveton Circle
Sparks, MD 21152
Telephone.: (410) 771 - 0100 or (800) 638 – 8663
- **Information Department:** Technical Service
- **Emergency telephone number:**
In case of a chemical emergency, spill, fire, exposure, or accident contact BD Diagnostic Systems (410) 771-0100 or (800)-638-8663, or ChemTrec at (800) 424-9300.



2 Composition/information on ingredients

- **Chemical characterization:** *Mixture*
- **Description:** *Mixture consisting of the following components.*

· **Dangerous Components:**

66-81-9	cycloheximide	1.22%
---------	---------------	-------

3 Hazards identification

- **Classification of the substance or mixture**
- **Classification according to Directive 67/548/EEC or Directive 1999/45/EC**
 **T; Toxic**
R61-25: May cause harm to the unborn child. Toxic if swallowed.
-  **Xn; Harmful**
R68: Possible risk of irreversible effects.
- **Classification system:**
The classification was made according to the latest editions of international substances lists, and expanded upon from company and literature data.
- **Label elements**
- **Labelling according to EU guidelines:**
The product has been classified and marked in accordance with regulations on hazardous materials.
- **Code letter and hazard designation of product:**
T Toxic
- **Hazard-determining components of labelling:**
cycloheximide
- **Risk phrases:**
61 May cause harm to the unborn child
25 Also toxic if swallowed.
68 Possible risk of irreversible effects.
- **Safety phrases:**
53 Avoid exposure - obtain special instructions before use.

(Contd. on page 2)

USA

Date Prepared: 03/01/2011

Reviewed On: 02/11/2011

Product Name: Dermatophyte Test Medium Base

(Contd. of page 1)

1/2 Keep locked up and out of the reach of children.

29/56 Do not empty into drains, dispose of this material and its container at hazardous or special waste collection point

36/37 Wear suitable protective clothing and gloves.

45 In case of accident or if you feel unwell, seek medical advice immediately.

60 This material and its container must be disposed of as hazardous waste.

· NFPA ratings (scale 0-4)


Health = 4

Fire = 0

Reactivity = 0

· HMIS ratings (scale 0-4)


Health = 4

Fire = 0

Reactivity = 0

4 First aid measures
· General information

Remove breathing apparatus only after contaminated clothing have been completely removed.

· **After skin contact** Immediately wash with water and soap and rinse thoroughly.

· **After eye contact** Rinse opened eye for 15 minutes under running water. Then consult a doctor.

· **After swallowing** Call a doctor immediately.

· **Information for doctor** Show this product label or this MSDS.

5 Firefighting measures
· Suitable extinguishing agents

CO₂, ABC multipurpose dry chemical or water spray. Fight larger fires with water spray or alcohol resistant foam.

· **Protective equipment:** No special measures required.

6 Accidental release measures

· **Personal precautions, protective equipment and emergency procedures** Not required.

· **Environmental precautions:** Wipe up with damp sponge or mop.

· **Methods and material for containment and cleaning up:** Ensure adequate ventilation.

· Reference to other sections

See Section 7 for information on safe handling.

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

7 Handling and storage
· Handling

· **Precautions for safe handling** Ensure good ventilation/exhaustion at the workplace.

· **Information about protection against explosions and fires:** No special measures required.

(Contd. on page 3)

Date Prepared: 03/01/2011

Reviewed On: 02/11/2011

Product Name: Dermatophyte Test Medium Base

(Contd. of page 2)

- **Storage**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Store away from oxidizing agents.
- **Further information about storage conditions:** Store in cool, dry conditions in well sealed containers.

8 Exposure controls/personal protection

- **Additional information about design of technical systems:** No further data; see Section 7.
- **Components with limit values that require monitoring at the workplace:**
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
- **Additional information:** The lists that were valid during the creation were used as basis.
- **Personal Protective Equipment**
- **General protective and hygienic measures**
Keep away from foodstuffs, beverages and feed.
Immediately remove all soiled and contaminated clothing.
Wash hands before breaks and at the end of work.
- **Breathing equipment:**
In case of brief exposure, use a chemical fume hood or a NIOSH/MSHA-approved respirator.
- **Protection of hands:**



Chemical resistant gloves (i.e. nitrile, or equivalent).

- **Eye protection:** Safety glasses
- **Body protection:** Protective work clothing (lab coat).

9 Physical and chemical properties

- | | |
|---------------------------------------------|--------------------------------------------------------------------------------------------|
| · General Information | |
| · Appearance: | |
| Form: | Solid. |
| Color: | Light beige |
| · Odor: | Characteristic |
| · Change in condition | |
| Melting point/Melting range: Not determined | |
| Boiling point/Boiling range: Not determined | |
| · Flash point: | Not applicable |
| · Flammability (solid, gaseous) | Product is not flammable. |
| · Danger of explosion: | Product is not explosive. However, formation of explosive air/vapor mixtures are possible. |
| · Density: | Not determined |
| · Solubility in / Miscibility with | |
| Water: | Soluble |
| · Solvent content: | |
| Organic solvents: | 0.0 % |

(Contd. on page 4)

Date Prepared: 03/01/2011

Reviewed On: 02/11/2011

Product Name: Dermatophyte Test Medium Base

(Contd. of page 3)

Solids content:	100.0 %
------------------------	---------

10 Stability and reactivity

- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Incompatible materials:** Incompatible material: strong oxidizers.
- **Hazardous decomposition products:** Nitrogen oxides (NO_x)

11 Toxicological information

- **Acute toxicity:**

- **LD/LC50 values that are relevant for classification:**

66-81-9 cycloheximide

Oral	LD50	2 mg/kg (rat)
------	------	---------------

- **Primary irritant effect:**

- **on the skin:** No irritating effect.

- **on the eye:** No irritating effect.

- **Sensitization:** No sensitizing effects known.

- **Subacute to chronic toxicity:**

Target organs: (Reproductive) central nervous system, eye/ear, musculoskeletal system, urogenital system.
(Carcinogenesis) Blood Leukemia, lymphomas, aplastic anemia.

- **Additional toxicological information:**

The product shows the following dangers according to internally approved calculation methods for preparations:
Toxic

12 Ecological information

- **Acquatic toxicity:** No further relevant information available.

- **Ecotoxicological effects:**

- **Other information:**

The ecological effects have not been thoroughly investigated, but currently none have been identified.

- **Additional ecological information:**

- **General notes:** Water hazard class 2 (Self-assessment): hazardous for water.

13 Disposal considerations

- **Waste treatment methods**

- **Recommendation**

Must not be disposed of with solid waste.

Dispose of material in accordance with federal (40 CFR 261.3), state and local requirements.

This product is not considered a RCRA hazardous waste.

- **Uncleaned packagings:**

- **Recommendation:** Disposal must be made according to state and federal regulations.

(Contd. on page 5)

Date Prepared: 03/01/2011

Reviewed On: 02/11/2011

Product Name: Dermatophyte Test Medium Base

(Contd. of page 4)

· **Recommended cleansing agent:** Water, if necessary with cleansing agents.

14 Transport information

· **DOT regulations:**

· **Hazard class:** -

· **Land transport ADR/RID (cross-border)**

· **ADR/RID class:** -

· **Maritime transport IMDG:**

· **IMDG Class:** -

· **Marine pollutant:** No

· **Air transport ICAO-TI and IATA-DGR:**

· **ICAO/IATA Class:** -

15 Regulatory information

· **SARA Section 355 (extremely hazardous substances)**

66-81-9 cycloheximide

· **SARA Section 313 (specific toxic chemical listings)**

None of the ingredients is listed.

· **TSCA (Toxic Substances Control Act)**

9002-18-0 agar

50-99-7 dextrose

143-74-8 phenol red

· **California Proposition 65 - Chemicals known to cause cancer**

None of the ingredients is listed.

· **California Proposition 65 - Chemicals known to cause reproductive toxicity for females:**

None of the ingredients is listed.

· **California Proposition 65 - Chemicals known to cause reproductive toxicity for males:**

None of the ingredients is listed.

· **California Proposition 65 - Chemicals known to cause developmental toxicity:**

66-81-9 cycloheximide

· **Carcinogenic categories**

· **NTP (National Toxicology Program)**

None of the ingredients is listed.

· **TLV (Threshold Limit Value established by ACGIH)**

None of the ingredients is listed.

· **National regulations**

· **Water hazard class:** Water hazard class 2 (Self-assessment): hazardous for water.

USA

(Contd. on page 6)



Material Safety Data Sheet
acc. to ISO/DIS 11014

Date Prepared: 03/01/2011

Reviewed On: 02/11/2011

Product Name: Dermatophyte Test Medium Base

(Contd. of page 5)

16 Other information

To the best of our knowledge, the information contained herein is accurate. However, neither Becton, Dickinson and Company or any of its subsidiaries assumes any liabilities whatsoever for the accuracy or completeness of the information contained herein. Final determination of suitability of any material is the sole responsibility of the user. All materials may present unknown hazards and should be used with caution. Although certain hazards are described herein, we can not guarantee that these are the only hazards that exist.

- **Department issuing MSDS:**
Environmental, Health & Safety
Created by Michael J. Spinazzola
- **Contact:** Technical Service Representative

USA