

## Certificate Of Analysis

---

### Product Information

---

**Product Number:** 27532501

**Product Name:** illustra™ MicroSpin™ G-25 Columns

**Pack Size:** 50 purifications

**Lot Number:** 9504943

**Expiry Date:** 20 July 2015

**Warning:** For research use only. Not recommended or intended for diagnosis of disease in humans or animals. Do not use internally or externally in humans or animals.

### Contents

---

28917922 illustra™ MicroSpin™ G-25 columns, Lot 9504943:  
1 × 50 pack

The columns contain Sephadex™ G-25 DNA Grade F and are supplied pre-equilibrated in double-distilled water containing 0.05% Kathon™ CG/ICP Biocide as a preservative.

28961634 Collection tubes  
1 × 50 pack

28963675 Caps  
1 × 50 pack

Nuclease testing: DNase/RNase tested.

Nuclease test results: No observable DNase/RNase.

### Conformance & Quality systems statement

---

This is to certify that this product conforms to GE Healthcare Life Sciences specifications.

All products are manufactured via a Quality System certified to **BS EN ISO9001:2008** and tested in accordance with documented quality procedures and approved as a result of meeting the required specification

### Electronic signature

---

This document has been electronically produced and is valid without a signature.

Version AM

[www.gelifesciences.com](http://www.gelifesciences.com)